

# Eyezen™ + 0, 1, 2 and 3 Eyefinity®

## ORDERING INSTRUCTIONS:

Internal use only - Do not distribute

Better Sight.  
Better Life.



**STEP 1**

**1 Select Vision Type:**  
Single Vision

**Lens Finishing** Lab Finishing

**Vision Type** Single Vision

**Material**

**Lens** Single Vision

**Lens Not Listed** Bifocal

**One Lens** Trifocal

**Balance Lens** Double  
Progressive  
Lenticular Single Vision

**Bevel** Lenticular Bifocal

**Bevel Not Listed** Lenticular Trifocal  
Quad  
Near Variable Focus

**Tint Color**

**Tint Not Listed**

**Tint Type**

**Base Curve**  
R  
L

**Edge**

**Edge Not Listed**

No Sample  Match Sample  
 Lighten  Darken  
Lighten or Darken by %

**STEP 2**

**2 Select desired Material.**  
Utilize the *Eyezen+* availability chart for full product availability.

**Lens Finishing** Lab Finishing

**Vision Type** Single Vision

**Material**

**Lens** Glass

**Lens Not Listed** Glass Hi Index

**One Lens** Plastic

**Balance Lens** Plastic Mid Index

**Bevel** Plastic Hi Index

**Bevel Not Listed** Polycarbonate  
Trivex

**Base Curve**  
R  
L

**Edge**

**Edge Not Listed**

**STEP 3**

**3** Scroll down in the Lens drop-down menu **below** “*Lab Choice*” selections and choose the desired *Eyezen+* product.

Remember, *Eyezen+* 0, 1, 2, and 3 come with different amounts of accommodative relief and are recommended for different age groups.

**Lens Finishing** Lab Finishing

**Vision Type** Single Vision

**Material** Plastic Hi Index

**Lens** Aspheric 1.67 w/ Crizal Avance UV AR (D) - Transitions Signature Gray

**Lens Not Listed** Aspheric 1.67 w/ Crizal Avance UV AR (D+UV) - Clear

**One Lens** Aspheric 1.67 w/ UNITY Classic AR (B) - Clear

**Balance Lens** Aspheric 1.67 w/ UNITY Elite AR (D) - Clear  
Aspheric 1.67 w/ UNITY Plus AR (C) - Clear

**Bevel** Aspheric 1.74 w/ Premium AR (C) - Clear

**Bevel Not Listed** Eyezen+ 1 1.67 - Clear  
Eyezen+ 1 1.67 - Polarized Brown C  
Eyezen+ 1 1.67 - Polarized Gray C  
Eyezen+ 1 1.67 - Polarized Green  
Eyezen+ 1 1.67 - Transitions Signature Brown  
Eyezen+ 1 1.67 - Transitions Signature Graphite Green  
Eyezen+ 1 1.67 - Transitions Signature Gray

**Base Curve**  
R  
L

**Sample ken**  
%

**STEP 4**

**4** Always include an Optical Center Height in your measurements and complete the order as usual.\*

Visit [EyezenPro.com](http://EyezenPro.com) to learn more.

\*Referred to as fitting height or height in VisionWeb.



## EYEZEN™+ LENSES AS A FULLY-COVERED BENEFIT!

Did you know VSP Vision Care patients fully covered for polycarbonate can choose Eyezen+ polycarbonate lenses as their fully covered option!

Make sure you talk to your eligible VSP patients about this great opportunity!

## DISPENSING INFORMATION

There are four unique Eyezen+ lens designs, and the patient's age will determine which Eyezen+ design is right for them. The difference between the designs is the amount of accommodative relief the patient will receive in the lens:



**EYEZEN+ 0**

**Patients age 17 and younger**  
(0.00 diopters)



**EYEZEN+ 1**

**Patients age 18 to 34**  
(0.40 diopters)



**EYEZEN+ 2**

**Patients age 35 to 44**  
(0.60 diopters)



**EYEZEN+ 3**

**Patients age 45 to 50\***  
(0.85 diopters)

UNDER 18



AGE 50

*Patient's age will determine which design is right for them.*

Eyezen+ lenses are fitted, ordered, and dispensed like a traditional single vision lens, except Eyezen+ lenses also **require a fitting height**. Simply measure from the center of the pupil to the lowest point of the lens to acquire fitting height.

\*It is recommended you consider progressive lenses with your patients age 45-50.

## RECOMMEND EYEZEN+ LENSES TODAY!

Go to [EyezenPro.com](http://EyezenPro.com) for additional information and materials

\*For Eyezen+ designs 1,2, and 3 containing accommodative relief.

\*\*Eyezen+ lenses block at least 20% of Harmful Blue Light, which is the high energy waves found between 415-455nm (blue-violet light).

©2017 Essilor of America, Inc. All Rights Reserved. Unless indicated otherwise, all trademarks are the property of Essilor International.  
LEZN0000XX PDF/HB 9/17

Eyezen™+